



About

The SWISSCI LCP Modular range aims to offer more flexibility to the membrane protein crystallographers LCP experiment. There are three choices of LCP Bases which are combined with one of five pre-affixed LCP Double Sided Pressure Adhesives. The Pressure Adhesives are of varying thickness and well dimensions. To complete the LCP sandwich plate there is a choice between two sizes of Cover Glass or three UVP Cover Films that are of varying thickness.



FACTS

- 5 unique LCP adhesives
- 3 bases
- 5 covers
- 75 different combinations

FEATURES

- Multiple configurations
- UV transparent bases
- Wider range of experimental options

LCP Sandwich Options and Codes

LCP-SBSG-60	Swissci Modular SBS XPOL Glass Base + 60µm 5.5mm Double Sided Adhesive
LCP-SBSG-60.7	Swissci Modular SBS XPOL Glass Base + 60µm 7mm Double Sided Adhesive
LCP-SBSG-100	Swissci Modular SBS XPOL Glass Base + 100µm 5.5mm Double Sided Adhesive
LCP-SBSG-120.7	Swissci Modular SBS XPOL Glass Base + 120µm 7mm Double Sided Adhesive
LCP-SBSG-200	Swissci Modular SBS XPOL Glass Base + 200µm 5.5mm Double Sided Adhesive
LCP-BASEG-60	Swissci Modular Glass Base + 60µm 5.5mm Double Sided Adhesive
LCP-BASEG-60.7	Swissci Modular Glass Base + 60µm 7mm Double Sided Adhesive
LCP-BASEG-100	Swissci Modular Glass Base + 100µm 5.5mm Double Sided Adhesive
LCP-BASEG-120.7	Swissci Modular Glass Base + 120µm 7mm Double Sided Adhesive
LCP-BASEG-200	Swissci Modular Glass Base + 200µm 5.5mm Double Sided Adhesive
LCP-UVXPO-60	Swissci Modular UVXPO Base + 60µm 5.5mm Double Sided Adhesive
LCP-UVXPO-60.7	Swissci Modular UVXPO Base + 60µm 7mm Double Sided Adhesive
LCP-UVXPO-100	Swissci Modular UVXPO Base + 100µm 5.5mm Double Sided Adhesive
LCP-UVXPO-120.7	Swissci Modular UVXPO Base + 120µm 7mm Double Sided Adhesive
LCP-UVXPO-200	Swissci Modular UVXPO Base + 200µm 5.5mm Double Sided Adhesive
LCP-G175	Swissci Modular Cover Glass 73x109mm NO 1.5
LCP-G175L	Swissci Modular Cover Glass 80x115mm NO 1.5
LCP-UV40	Swissci Modular UVP Cover Film 40 Micron 109x73mm
LCP-UV100	Swissci Modular UVP Cover Film 100 Micron 109x73mm
LCP-UV200	Swissci Modular UVP Cover Film 200 Micron 109x73mm
LCP-TRAYW	Swissci Modular 96 Well Tray Moulding Plate White PS
LCP-TRAYW-LT	Swissci Modular 96 Well Tray Moulding Plate White PS Low Tack
LCP-PADDLE	Swissci Modular Swissci Sealing Paddle

Order information

Item:

LCP Modular Range

Codes:

See table

Details:

Pack of 10

Contact: sales@swissci.com

Web: www.swissci.com