



## About

The SWISSCI LCP Modular range aims to offer more flexibility to the membrane protein crystallographers LCP experiment. There are three choices of LCP Bases which are combined with one of five pre-affixed LCP Double Sided Pressure Adhesives. The Pressure Adhesives are of varying thickness and well dimensions. To complete the LCP sandwich plate there is a choice between two sizes of Cover Glass or three UVP Cover Films that are of varying thickness.



## FACTS

- 5 unique LCP adhesives
- 3 bases
- 5 covers

## FEATURES

- Multiple configurations
- UV Transparent bases
- Wider range of experimental options

## LCP Sandwich Options and Codes

**LCP Bases:** SBS XPOL Glass Base, Glass Base and UVXPO Plastic Base  
**Spacers:** 60 µm, 100 µm, 120 µm, 200 µm. **Well sizes:** 5.5 mm and 7 mm.  
**Covers:** Cover Glass, Large Cover Glass, 40 µm Cover Film, 100 µm Cover Film, 200 µm Cover film covers

LCP-SBSG-60	Swissci Modular SBS XPOL Glass Base + 60µm 5.5mm Double Sided Adhesive (pack 10)
LCP-SBSG-60.7	Swissci Modular SBS XPOL Glass Base + 60µm 7mm Double Sided Adhesive (pack 10)
LCP-SBSG-100	Swissci Modular SBS XPOL Glass Base + 100µm 5.5mm Double Sided Adhesive (pack 10)
LCP-SBSG-120.7	Swissci Modular SBS XPOL Glass Base + 120µm 7mm Double Sided Adhesive (pack 10)
LCP-SBSG-200	Swissci Modular SBS XPOL Glass Base + 200µm 5.5mm Double Sided Adhesive (pack 10)
LCP-BASEG-60	Swissci Modular Glass Base + 60µm 5.5mm Double Sided Adhesive (pack of 10)
LCP-BASEG-60.7	Swissci Modular Glass Base + 60µm 7mm Double Sided Adhesive (pack 10)
LCP-BASEG-100	Swissci Modular Glass Base + 100µm 5.5mm Double Sided Adhesive (pack 10)
LCP-BASEG-120.7	Swissci Modular Glass Base + 120µm 7mm Double Sided Adhesive (pack 10)
LCP-BASEG-200	Swissci Modular Glass Base + 200µm 5.5mm Double Sided Adhesive (pack 10)
LCP-UVXPO-60	Swissci Modular UVXPO Base + 60µm 5.5mm Double Sided Adhesive (pack of 10)
LCP-UVXPO-60.7	Swissci Modular UVXPO Base + 60µm 7mm Double Sided Adhesive (pack 10)
LCP-UVXPO-100	Swissci Modular UVXPO Base + 100µm 5.5mm Double Sided Adhesive (pack 10)
LCP-UVXPO-120.7	Swissci Modular UVXPO Base + 120µm 7mm Double Sided Adhesive (pack 10)
LCP-UVXPO-200	Swissci Modular UVXPO Base + 200µm 5.5mm Double Sided Adhesive (pack 10)
LCP-G175	Swissci Modular Cover Glass 73x109mm NO 1.5 (pack of 10)
LCP-G175L	Swissci Modular Cover Glass 80x115mm NO 1.5 (pack of 10)
LCP-UV40	Swissci Modular UVP Cover Film 40 Micron 109x73mm (pack of 10)
LCP-UV100	Swissci Modular UVP Cover Film 100 Micron 109x73mm (pack of 10)
LCP-UV200	Swissci Modular UVP Cover Film 200 Micron 109x73mm (pack of 10)
LCP-TRAYW	Swissci Modular 96 Well Tray Moulding Plate White PS (pack of 10)
LCP-TRAYW-LT	Swissci Modular 96 Well Tray Moulding Plate White PS Low Tack (pack of 10)
LCP-PADDLE	Swissci Modular Swissci Sealing Paddle (pack of 10)

**Contact:** sales@swissci.com

**Tel:** +44 (0)1932 481444

**Web:** www.swissci.com