

Wizard PEG Ion 48 salt screen with PEG 8K and 10K

(96 formulations, 1.7 ml each in a 96-well block plate)

1008656

Well	Precipitation Reagent	Salt (0.2M)
A1	20% (w/v) PEG 8000	Ammonium acetate
A2	20% (w/v) PEG 8000	Ammonium chloride
A3	20% (w/v) PEG 8000	Ammonium fluoride
A4	20% (w/v) PEG 8000	Ammonium formate
A5	20% (w/v) PEG 8000	Ammonium iodide
A6	20% (w/v) PEG 8000	Ammonium nitrate
A7	20% (w/v) PEG 8000	Ammonium phosphate (monobasic)
A8	20% (w/v) PEG 8000	Ammonium phosphate (dibasic)
A9	20% (w/v) PEG 8000	Ammonium sulfate
A10	20% (w/v) PEG 8000	Ammonium sulfite
A11	20% (w/v) PEG 8000	Calcium acetate
A12	20% (w/v) PEG 8000	Calcium chloride
B1	20% (w/v) PEG 8000	Diammonium tartrate
B2	20% (w/v) PEG 8000	Formic acid potassium salt
B3	20% (w/v) PEG 8000	Lithium acetate
B4	20% (w/v) PEG 8000	Lithium chloride
B5	20% (w/v) PEG 8000	Lithium citrate
B6	20% (w/v) PEG 8000	Lithium nitrate
B7	20% (w/v) PEG 8000	Lithium sulfate
B8	20% (w/v) PEG 8000	Magnesium acetate
B9	20% (w/v) PEG 8000	Magnesium chloride
B10	20% (w/v) PEG 8000	Magnesium formate
B11	20% (w/v) PEG 8000	Magnesium nitrate
B12	20% (w/v) PEG 8000	Magnesium sulfate
C1	20% (w/v) PEG 8000	Potassium acetate
C2	20% (w/v) PEG 8000	Potassium chloride
C3	20% (w/v) PEG 8000	Potassium fluoride
C4	20% (w/v) PEG 8000	Potassium iodide
C5	20% (w/v) PEG 8000	Potassium nitrate
C6	20% (w/v) PEG 8000	50mM Potassium phosphate (dibasic)
C7	20% (w/v) PEG 8000	Potassium phosphate (monobasic)
C8	20% (w/v) PEG 8000	Potassium sulfate
C9	20% (w/v) PEG 8000	Potassium thiocyanate
C10	20% (w/v) PEG 8000	Potassium/sodium tartrate
C11	20% (w/v) PEG 8000	Sodium acetate
C12	20% (w/v) PEG 8000	150mM Sodium chloride
D1	20% (w/v) PEG 8000	Sodium citrate
D2	20% (w/v) PEG 8000	Sodium fluoride
D3	20% (w/v) PEG 8000	Sodium formate
D4	20% (w/v) PEG 8000	Sodium iodide
D5	20% (w/v) PEG 8000	Sodium isothiocyanate
D6	20% (w/v) PEG 8000	Sodium nitrate
D7	20% (w/v) PEG 8000	Sodium phosphate (dibasic)
D8	20% (w/v) PEG 8000	Sodium phosphate (monobasic)
D9	20% (w/v) PEG 8000	Sodium sulfate
D10	20% (w/v) PEG 8000	Sodium tartrate (dibasic)
D11	20% (w/v) PEG 8000	Tripotassium citrate
D12	20% (w/v) PEG 8000	Zinc acetate

Well	Precipitation Reagent	Salt (0.2M)
E1	10% (w/v) PEG 10000	Ammonium acetate
E2	10% (w/v) PEG 10000	Ammonium chloride
E3	10% (w/v) PEG 10000	Ammonium fluoride
E4	10% (w/v) PEG 10000	Ammonium formate
E5	10% (w/v) PEG 10000	Ammonium iodide
E6	10% (w/v) PEG 10000	Ammonium nitrate
E7	10% (w/v) PEG 10000	Ammonium phosphate (monobasic)
E8	10% (w/v) PEG 10000	Ammonium phosphate (dibasic)
E9	10% (w/v) PEG 10000	Ammonium sulfate
E10	10% (w/v) PEG 10000	Ammonium sulfite
E11	10% (w/v) PEG 10000	Calcium acetate
E12	10% (w/v) PEG 10000	50mM Calcium chloride
F1	10% (w/v) PEG 10000	Ammonium tartrate dibasic
F2	10% (w/v) PEG 10000	Formic acid potassium salt
F3	10% (w/v) PEG 10000	Lithium acetate
F4	10% (w/v) PEG 10000	Lithium chloride
F5	10% (w/v) PEG 10000	Lithium citrate
F6	10% (w/v) PEG 10000	Lithium nitrate
F7	10% (w/v) PEG 10000	Lithium sulfate
F8	10% (w/v) PEG 10000	Magnesium acetate
F9	10% (w/v) PEG 10000	50mM Magnesium chloride
F10	10% (w/v) PEG 10000	Magnesium formate
F11	10% (w/v) PEG 10000	Magnesium nitrate
F12	10% (w/v) PEG 10000	Magnesium sulfate
G1	10% (w/v) PEG 10000	Potassium acetate
G2	10% (w/v) PEG 10000	100mM Potassium chloride
G3	10% (w/v) PEG 10000	Potassium fluoride
G4	10% (w/v) PEG 10000	Potassium iodide
G5	10% (w/v) PEG 10000	Potassium nitrate
G6	10% (w/v) PEG 10000	Potassium phosphate (dibasic)
G7	10% (w/v) PEG 10000	100mM Potassium phosphate (monobasic)
G8	10% (w/v) PEG 10000	100mM Potassium sulfate
G9	10% (w/v) PEG 10000	Potassium thiocyanate
G10	10% (w/v) PEG 10000	Potassium/Sodium tartrate
G11	10% (w/v) PEG 10000	100mM Sodium acetate
G12	10% (w/v) PEG 10000	Sodium chloride
H1	10% (w/v) PEG 10000	50mM Sodium citrate
H2	10% (w/v) PEG 10000	Sodium fluoride
H3	10% (w/v) PEG 10000	Sodium formate
H4	10% (w/v) PEG 10000	Sodium iodide
H5	10% (w/v) PEG 10000	Sodium isothiocyanate
H6	10% (w/v) PEG 10000	Sodium nitrate
H7	10% (w/v) PEG 10000	Sodium phosphate (dibasic)
H8	10% (w/v) PEG 10000	50mM Sodium phosphate (monobasic)
H9	10% (w/v) PEG 10000	Sodium sulfate
H10	10% (w/v) PEG 10000	Sodium tartrate (dibasic)
H11	10% (w/v) PEG 10000	100mM Potassium citrate tribasic
H12	10% (w/v) PEG 10000	Zinc acetate