

## About

The Diatube™ is a intuitive all-in-one device for dialysis of up to 20mls of sample solution. The Diatube™ requires no preparation; it is ready to use and is ideal for large sample dialysis.

The membrane is made from regenerated cellulose and is specially glued to the device to form a leak free seal; no knotting, folding or closures required.

The Swissci Diatube™ has a MWCO of 10kDa.



### FACTS

- Regenerated cellulose membrane
- 20 ml capacity
- Multiple MWCO's available

### FEATURES

- All-in-one device
- Ready to use
- Rapid Dialysis

## Product Specifications

- Ready to use, just add sample
- Ideal for large sample volumes up to 20 ml
- No cross contamination or leakage
- Regenerated cellulose membrane; sulphur and heavy metal free
- Molecular weight cut off of 10,000 Daltons
- Excellent sample recovery without disturbing valuable protein complexes
- Large surface area for rapid dialysis
- Gentle and efficient desalting
- Caps to prevent spillage/evaporation
- Single use - disposable device; no cleaning required

### Order Information

**Item:** Diatube 20 ml Dialyser

**Code:** W92003; W92005;

W92010; W92030

**Details:** Pack of 10

**Contact:** sales@swissci.com

**Web:** www.swissci.com