


About

Diacon™ is a low volume dialyser for protein solutions of up to 0.5ml.

The Swissci Diacon™ uses specially selected regenerated cellulose membranes for dialysis in a provided floating rack for 2-3 hours.



Product Specifications

- Ready to use, just add sample
- Ideal for small sample volumes up to 0.5ml
- No cross contamination or leakage
- Regenerated cellulose membrane; sulphur and heavy metal free
- Molecular weight cut offs of 3, 5, 10 or 30 KDa available
- Excellent sample recovery without disturbing valuable protein complexes
- Dialysis complete in 2-3 hours
- Gentle and efficient desalting
- Caps to prevent spillage/evaporation
- Floating tube rack provided holds up to 18 Diacon at any one time
- Single use - disposable device; no cleaning required

FACTS

- Regenerated cellulose membrane
- 0.5 ml capacity
- Multiple MWCO's available

FEATURES

- Dialysis in 2-3 hours
- Floating rack included
- Ready to use

Order Information

Item: Diacon 0.5 ml Dialyser

Codes: W72003; W72005;
W72010; W72030.

Details: Pack of 40

Contact: sales@swissci.com

Web: www.swissci.com