



About

The SWISSCI LCP Modular range aims to offer more flexibility to the membrane protein crystallographers LCP experiment. There are three choices of LCP Bases which are combined with one of five pre-affixed LCP Double Sided Pressure Adhesives. The Pressure Adhesives are of varying thickness and well dimensions. To complete the LCP sandwich plate there is a choice between three sizes of Cover Glass or three UVP Cover Films that are of varying thickness.



FEATURES

- 5 unique LCP Adhesives
- 3 Bases
- 6 Covers
- 90 different experimental options
- UV transparent Bases + Covers
- Hydrophobic option
- All options SBS compatible
- Bespoke film cutting possible
- Bespoke Adhesive well design possible

LCP Sandwich Options and Codes

LCP-SBSG-60	Swissci Modular SBS XPOL Glass Base + 60μm 5.5mm Double Sided Adhesive
LCP-SBSG-60.7	Swissci Modular SBS XPOL Glass Base + 60μm 7mm Double Sided Adhesive
LCP-SBSG-100	Swissci Modular SBS XPOL Glass Base + 100μm 5.5mm Double Sided Adhesive
LCP-SBSG-120.7	Swissci Modular SBS XPOL Glass Base + 120μm 7mm Double Sided Adhesive
LCP-SBSG-200	Swissci Modular SBS XPOL Glass Base + 200μm 5.5mm Double Sided Adhesive
LCP-BASEG-60	Swissci Modular Glass Base + 60μm 5.5mm Double Sided Adhesive
LCP-BASEG-60.7	Swissci Modular Glass Base + 60μm 7mm Double Sided Adhesive
LCP-BASEG-100	Swissci Modular Glass Base + 100μm 5.5mm Double Sided Adhesive
LCP-BASEG-120.7	Swissci Modular Glass Base + 120μm 7mm Double Sided Adhesive
LCP-BASEG-200	Swissci Modular Glass Base + 200μm 5.5mm Double Sided Adhesive
LCP-UVXPO-60	Swissci Modular UVXPO Base + 60μm 5.5mm Double Sided Adhesive
LCP-UVXPO-60.7	Swissci Modular UVXPO Base + 60μm 7mm Double Sided Adhesive
LCP-UVXPO-100	Swissci Modular UVXPO Base + 100μm 5.5mm Double Sided Adhesive
LCP-UVXPO-120.7	Swissci Modular UVXPO Base + 120μm 7mm Double Sided Adhesive
LCP-UVXPO-200	Swissci Modular UVXPO Base + 200μm 5.5mm Double Sided Adhesive
LCP-G100	Swissci Modular Cover Glass 73x109mm 100 μm
LCP-G175	Swissci Modular Cover Glass 73x109mm 175 μm
LCP-G175L	Swissci Modular Cover Glass 80x115mm 175 μm
LCP-UV40	Swissci Modular UVP Cover Film 109x73mm 40 μm
LCP-UV100	Swissci Modular UVP Cover Film 109x73mm 100 μm
LCP-UV200	Swissci Modular UVP Cover Film 109x73mm 200 μm
LCP-TRAY	Swissci Modular 96 Well SBS Tray holder
LCP-PADDLE	Swissci Modular Swissci Sealing Paddle

NOTE: All Bases and Adhesives are also available separately. Message sales@swissci.com to find out more.

Order information

Item:

LCP Modular Range

Codes:

See table

Details:

Packs of 10

Contact: sales@swissci.com

Web: www.swissci.com