


## About

The SWISSCI LCP Modular range aims to offer more flexibility to the membrane protein crystallographers LCP experiment. There are three choices of LCP Bases which are combined with one of seven pre-affixed LCP Double Sided Pressure Adhesives.


The Pressure Adhesives are of varying thickness (60-800  $\mu\text{m}$ ) and well dimensions (5.5-7mm). To complete the LCP sandwich plate there is a choice between three sizes of Cover Glass varying in thickness (100-175  $\mu\text{m}$ ) or four UVP Cover Films (13-200  $\mu\text{m}$ ).



## LCP Sandwich Options and Codes

<b>LCP-SBSG-60/60.7/100/120.7/200/400/800</b>	E.g. Swissci Modular SBS XPOL Hydrophobic Glass Base + 60 $\mu\text{m}$ 5.5mm Double Sided Adhesive (pk 10)
<b>LCP-BASEG-60/60.7/100/120.7/200/400/800</b>	E.g. Swissci Modular Glass Base + 60 $\mu\text{m}$ 5.5mm Double Sided Adhesive (pk 10)
<b>LCP-UVXPO-60/60.7/100/120.7/200/400/800</b>	E.g. Swissci Modular UVXPO Base + 60 $\mu\text{m}$ 5.5mm Double Sided Adhesive (pk 10)
<b>LCP-G100</b>	Swissci Modular Cover Glass 73x109mm 100 $\mu\text{m}$ (pk 10)
<b>LCP-G175</b>	Swissci Modular Cover Glass 73x109mm 175 $\mu\text{m}$ (pk 10)
<b>LCP-G175L</b>	Swissci Modular Cover Glass 80x115mm 175 $\mu\text{m}$ (pk 10)
<b>LCP-UV13</b>	Swissci Modular UVP Cover Film 109x73mm 13 $\mu\text{m}$ (pk 10)
<b>LCP-UV40</b>	Swissci Modular UVP Cover Film 109x73mm 40 $\mu\text{m}$ (pk 10)
<b>LCP-UV100</b>	Swissci Modular UVP Cover Film 109x73mm 100 $\mu\text{m}$ (pk 10)
<b>LCP-UV200</b>	Swissci Modular UVP Cover Film 109x73mm 200 $\mu\text{m}$ (pk 10)
<b>LCP-TRAY</b>	Swissci Modular 96 Well SBS Tray holder (pk 10)
<b>LCP-PADDLE</b>	Swissci Modular Swissci Sealing Paddle (pk 5)

NOTE: All Bases and Adhesives are also available separately. Should you require a different adhesive thickness please let us know as we may be able to accommodate your requirements. Please message [sales@swissci.com](mailto:sales@swissci.com) for more information.



## FEATURES

- 7 unique LCP Adhesives
- 3 Bases
- 7 Covers
- Over 147 experimental options
- UV transparent Bases + Covers
- Hydrophobic option
- All options SBS compatible
- Bespoke film cutting possible
- Bespoke Adhesive well design possible

## Order information

### Item:

LCP Moldular Range

### Codes:

See table

### Details:

Packs of 10

**Contact:** [sales@swissci.com](mailto:sales@swissci.com)

**Web:** [www.swissci.com](http://www.swissci.com)